Gebruiksaanweisung

Multisil-Epithetik Fasern

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Instructions for use

Multisil-Epithetik fibres

Before using the product, please read these instructions for careful use.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Multisil-Epithetik fibres are viscose fibres available in silver, white, beige, signal red, ruby coloured, purple, bordeaux, blue, mocha and ochre.

2. Indication range

Multisil-Epithetik fibres are used to reproduce blood vessels or other characteristic skin features by mixing them with Multisil-Epithetik transparent or Multisil-Epithetik city/country/beach as well as to recreate downy hair, thin veins on the surface or similar by applying them onto fresh Multisil sealing agent. Multisil-Epithetik fibres provide the epithesis with a natural appearance.

3. Contraindication

Based on the current state of knowledge there are no known contraindications.

4. Safety recommendations and hazard warning

In case of skin contact, rinse with copious amount of water. Please read the material safety datasheets for the product (request separately).

5. Storage and durability

5.1 Storage

The fibres must be stored at room temperature (not above 28° C) under dry conditions; the container (bottle) must be sealed well after the use.

5.2 Durability

The fibres feature unlimited durability if they are stored in the original, sealed container at room temperature (not above 28° C).

6. Processing

As required, the Multisil-Epithetik fibres can be mixed into the epithesis material or spread onto the sealing agent that is still wet.

7. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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GmbH & Co. KG · Weissenhorner Str. 2 · 89250 Senden · Germany